

Pocket Peak®

Mechanical Peak Flow Meter



Redefining Patient Care Beyond Expectations

Small, lightweight, and easy to use, Pocket Peak® is designed to deliver results you can trust.

Redefining Patient Care — Diagnostic Confidence

Manufactured to tight tolerances with rugged polymer construction, Pocket Peak is Best-in-Class for peak expiratory flow measurement.

Redefining Patient Care — Improving Outcomes

As a part of our family of asthma and disease management products, Pocket Peak embodies years of dedication to advancing the quality of lung health monitoring at home.

We continue redefining patient care so that you can achieve improved patient outcomes, more efficiently, more profitably.


Redefining Patient Care — Proven Technology

- **Zone Monitoring System.** RYG Color Zones. Easily adjusted for personal settings.
- **Easy to Clean.** No disassembly required.
- **Patented Design.** Precision engineered from the highest quality materials providing a lightweight, durable, pocket-sized design.
- **Universal Mouthpiece.** Sized for both adults and children and can be stowed conveniently on the side of the meter.
- **Compliant.** Each meter is ATS/EU scale compliant and is individually calibrated to ensure accurate patient test results.



Redefining Accuracy. Beyond Expectations.

Each Pocket Peak is individually calibrated, bringing you and your patient the most accurate peak flow measurements economically.



Technical Specifications

Range:

50-720 l/min BTPS

Scale graduation increment:

10 l/min

Accuracy:

±10% of reading

Repeatability:

within 10 l/min

Inter-Device Variability:

within 5%

Dimensions:

L,W,D: 86 x 81 x 30 mm

Material:

DOW ABS Magnum 9010 or Delta Polymers

ABS 650

Weight:

45 grams (inc. mouthpiece)

Color:

White

Performance Compliance:

ATS 1994 Standard, 1997 NAEPP (National Asthma Educator Prevention Program)

Includes:

Multilingual instructions; 60 page Asthma Diary Chart; Predicted set (Adult/Pediatric)



©Copyright nSpire Health 2007. Due to continual innovations, nSpire Health reserves the right to change specifications without notice. Pocket Peak is a registered trademark of nSpire Health.

Contact Information

Phone: 800.574.7374 303.666.5555

Email: sales@nspirehealth.com

Web: www.nspirehealth.com

PN 633078 REV A