


Disposable Surgical Pencil


- Double sterilization pouches
 - Peeling pouches
 - Sterilization marker on the pouch
 - Traceability label for the patient's file
 - Quality control : 100%
- Overmolded connectors
 - Extra flexible cables
 - Anti-rotation system blade



<i>Description</i>	<i>Part Number</i>	<i>Connector</i>	<i>Cable length</i>
Hand control - Blade electrode	BE 100	Valley lab	3 m
Foot control - Blade electrode	BE 500	Valley lab	3 m

DISPOSABLE ELECTRODES

<i>Part Numbers</i>	<i>Description</i>	<i>Length</i>
MI LA70	Blade electrode	70 mm
MI LA150	Blade electrode	150 mm
MI A170	Needle electrode	70 mm
MI AI150	Needle electrode	150 mm
MI BO70	Ball electrode	70 mm



Disposable Electrosurgical Plate



PLACOM[®]

Round Electrosurgical plate

Innovation

- Waterproofed, flexible and extra-thin Woven backing, the plate adapts itself to the skin.
- The thin conductive part improves Again the plate conformability.
- The round shape eases the plate positioning on the patient and enables different orientations.
- Optimized conductive surface for an excellent current dispersion in order to obtain a very low temperature (53°C)
- Size for adult and child (from 7 kg)
- Traceability sticker for a secured follow-up
- Compatible with the main HF generators
- The plate is hypoallergenic and does not contain latex

Technical characteristics :

- Adhesive surface : 167 cm²
- Conductive surface : 107.5 cm²
- Maximum impedance : $\geq 3 \Omega$ (according to AAMI standards)
- Maximum temperature : 30C
- Conservation : 24 months

Part numbers	Qty / pouch	Qty / Box
MF 3.05.5000	1	100
MF 3.05.5005	5	150

Connection cables - Length : 4 m



REM Valleylab
MI 3.05. VALD.40



Jack 6,35 standard
MI 3.05.BLUE.40



Jack 6,35 Erbe
MI 3.05.BLRB.40

Distributed by :



JK Medical Systems Pvt. Ltd.
1052 Munusamy Salai, West K.K. Nagar,
Chennai - 600 078, India.
Tel : 044 2474 2615, 2371 0009
Fax : +91 44 2474 4419
email : jkmed@vsnl.com,
Web : www.jkmedical.co.in

comepa
industries

88, Avenue du General de Gaulle - 93541 - Bagnolet Cedex
Tel. : +33 (0)1 42 43 50 50 - Fax : +33 (0)1 42 43 64 44
E-mail : info@comepa.com - www.comepa.com

Surgery