

Reusable Surgical Pencil

HAND CONTROL REUSABLE SURGICAL PENCIL

Autoclavable



Autoclavable at 134°C during 18 minutes
Change the pencil every 60 cycles

Safety

Reinforced cable at pencil & connector's junction

Reliability


Hexagonal anti-rotation system
100% single tested
Waterproofed switches




FINGER SWITCH ADAPTORS



Berchtold / Martin generator



Erbe ACC/ICC generator with 5 mm pin




Erbe T series generator with 4 mm pin

<i>Part Number</i>	<i>Description</i>	<i>Cable</i>	<i>Pack</i>
PMS 950	Autoclavable Surgical Pencil	5m	1
PMS 800	Finger Switch Adoptor for Berchtold/Martin generators	-	1
PMS 850	Finger Switch Adoptor for Erbe ACC/ICC with 5mm pin	-	1
PMS 860	Finger Switch Adoptor for Erbe T series with 4mm pin	-	1


Reusable Electrodes

REUSABLE BALL ELECTRODES




<i>Part Number</i>	<i>Description</i>	<i>Pack</i>
PMS1410/1403/1435	5cm Ball electrode 3mm/4mm/5mm	5
PMS1408/1407/1420	11cm Ball electrode 5mm/4mm/5mm	3
PMS1415/1416/1417	15cm Ball electrode 3mm/4mm/5mm	3

REUSABLE BLADE ELECTRODES



<i>Part Number</i>	<i>Description</i>	<i>Pack</i>
PMS1001	7 cm Blade Electrode	5
PMS1008	10 cm Blade Electrode	3
PMS1006	18 cm Blade Electrode	3

REUSABLE LOOP ELECTRODE



<i>Part Number</i>	<i>Description</i>	<i>Pack</i>
PMS1301/1302/1303/ 1304/1305/1306/1307	4.5cm Re-usable Loop electrodes 3mm/6mm/8mm/10mm/ 12mm/15mm/20mm	5
PMS1311/1312/1313/ 1314/1315/1316	10cm Re-usable Loop electrodes 6mm/8mm/10mm/ 12mm/15mm/20mm	3



Distributed by :

JK Medical Systems Pvt. Ltd.

1052 Munusamy Salai, West K.K. Nagar, Chennai - 600 078, India.

Tel : 044 2474 2615, 2371 0009 Fax : +91 44 2474 4419

email : jkmed@vsnl.com Web : www.jkmedical.co.in

prima
Medical Limited
www.prima-medical.com