

LungChek™

PEAK FLOW METER



HIGH STANDARD ACCURACY

The high standard accuracy peak flow measurement device for Asthma management, 'Lungchek™' enables objective measurement of airway constriction, suggesting actions necessary for Asthma management. Individually calibrated & tested for absolute accuracy.

INDICATION FOR USE

- Asthma patients to monitor airway constriction and lung function.
- Healthcare professionals assess respiratory conditions.
- Clinical settings for regular lung function monitoring.

TECHNICAL SPECIFICATIONS

Range	60-800 l/min (for adults) 60-600 l/min (for pead)
Scale Increment	10 l/min
Accuracy	< +/- 5% or 10 l/min
Repeatability	< +/- 3%
Calibration	100% (individual unit)
Weight	60gms
Dimension	Height - 150mm Diameter - 42mm

FEATURES

Large, high visibility scale
Available with European/ATS scale
Asthma management color zones
Customisation possible for promotional purposes
Meets and exceeds all international performance standards
Available colors: Blue & Orange

Cat No. 6051/Adult, 6052/Pead

Pkts. 01x20 Nos.



**JK MEDICAL
SYSTEMS PVT LTD**

1068, Munusamy Salai, West K K Nagar, Chennai - 600078, Tamilnadu, India.
Phone : +91 44 2366 4761 / 62, Email : sales@jkmedical.co.in
www.jkmedical.co.in