

INDSPIRE™

NEBULIZER KIT



COMFORT
RELIABILITY
&
EFFICIENCY

The INDSPIRE NEBULIZER KIT is designed for effective delivery of aerosolized medication, ensuring optimal therapeutic outcomes. Engineered for patient comfort, reliability, and efficiency, it is suitable for routine respiratory treatments at home or in clinical settings.

PERFORMANCE CHARACTERISTICS

- Efficient Medication Delivery: Fine aerosol particles (<5 µm) ensure deep lung penetration.
- Optimized Treatment Duration: Completes a full session in just 8–12 minutes.
- High Breathable Fraction: 63% of particles reach the lower respiratory tract for maximum therapeutic effect.
- Durable & Safe: Anti-crush tubing ensures safe handling and long-term use.

BENEFITS

- Compact, lightweight, and easy to use
- Suitable for both adults and children
- Quiet operation for a comfortable treatment experience
- Compatible with standard nebulizer medications

TECHNICAL SPECIFICATIONS

Chamber Volume	5 ml
Particle Size (MMD)	< 3.9 µm
Nebulization Output Rate (Minimum)	250 mg/min
Nebulization Output Rate (Maximum)	520 mg/min
Breathable Fraction < 5 µm (FPF)	63%
Tubing	Approximately 2m anti-crush, medical-grade PVC tubing
Typical Nebulization Duration	8 to 12 minutes per treatment session

ORDER INFORMATION

Cat No	Kit with Adult Mask	Kit with PAED Mask	Kit with Mouthpiece	Kit with Adult & PAED Mask	MOQ
2001	■				100 Nos
2002		■			100 Nos
2003			■		100 Nos
2004				■	50 Nos