


EASYFLO

DOUBLE WATER TRAP CIRCUIT (ADULT)



**JK MEDICAL
SYSTEMS PVT LTD**

Easyflo ventilator circuit is a medical device used to deliver oxygen, remove carbon dioxide from the patient. It comes with different configurations depending on customer's requirement.



Features

- Successful notch design eliminates the risk of disconnection and leakage
- Fully customizable as per user's demand
- All circuits are manufactured in class 7 standard clean room
- All circuits are comply with ISO13485:2016 Quality regulation and GMP standard
- Specially designed to use in ICU
- Sterile option available

Specifications

- 4 x 750 cm – 22 cm Corrugated Hose
- 1 x 600 cm – 22 cm Corrugated Hose
- 4 x 22F end Connectors
- 2 x water trap
- 1 x Y-Connector (with / without pressure ports)
- 1 x Elbow with cap

Order Information

Cat No	22mm/750 mm Tube	End Connectors	Parallel/ Bi-Wye Piece	Elbow with/without port	Water Trap	Catheter Mount	HME Filter	BV Filter	Neb kit	T-Piece	Proximal Line	M-M Connector	NIV Mask
1030	X	X	X	X	X								
1030/1	X	X	X	X	X	X							
1030/2	X	X	X	X	X		X						
1030/3	X	X	X	X	X		X	X					
1030/4	X	X	X	X	X	X	X						
1030/5	X	X	X	X	X	X	X	X					
1030/6	X	X	X	X	X				X	X	X	X	
1030/7	X	X	X	X	X					X			
1030/8	X	X	X	X	X							X	
1030/9	X	X	X	X	X		X	X					X



JK MEDICAL SYSTEMS PVT LTD

1068, Munusamy Salai, West K K Nagar, Chennai - 600078, Tamilnadu, India.

Phone : +91 44 2366 4761 / 62, Email : sales@jkmedical.co.in

www.jkmedical.co.in

